

Ultra Sensitive

High Tensile Strength

Superior Comfort







NEW GENERATION POWDER FREE NITRILE EXAMINATION GLOVES 2

Manufactured in accordance with Quality Management System ISO 13485 & 9001 Designed and Produced to meet the applicable requirements of EN455 (Parts 1, 2, 3 and 4) / EN374 (Parts 2 & 3)



Extra Sensitive



NEW GENERATION POWDER FREE NITRILE EXAMINATION GLOVES

Manufactured in accordance with Quality Management System ISO 13485 & 9001 Designed and Produced to meet the applicable requirements of EN455 (Parts 1, 2, 3 and 4) / EN374 (Parts 2 & 3)



APPROVED STANDARDS

These gloves have been tested and approved in accordance with the requirements of:

XXXXX

Medical Devices

European Standard EN455 Parts 1, 2, 3 and 4 European Directive 93/42/EEC



CE Class I Medical Devices

Personal Protective Equipment European Directive 89/686/EEC



C € PPE Cat III Complex Design

European Standard EN374 Parts 2 & 3







Food Handling European Directive EU No. 10/2011



STORAGE

Once opened contents should be shielded $+5^{\circ}C$ from direct sunlight, fluorescent lighting, x-rays, moisture and ozone.



DISPOSAL

Please ensure gloves are disposed of in an appropriate container. Do not reuse the gloves.

IMPORTANT

These gloves are 100% Nitrile (Acrylonitrile), and are designed to provide an excellent biological barrier and good resistance to a variety of chemicals, found in typical medical and industrial environments. For further information on chemical resistance, please contact our technical team on 01708 863868

This product does not contain natural rubber late: and is suitable for a latex free environment.

Made in Malaysia

Medicare Products Ltd, Unit B Dolphin Way, Purfleet RM19 1NZ, UK

AQL<1.5 (Level 2)

Performance Level	AQL	Inspection Levels
Level 3	< 0.65	G1
Level 2	<1.5	G1
Level 1	<4.0	S4

EN374-3:2003

Chemical	Level
Sodium Hydroxide 40%	6
Sulphuric Acid 50%	6
1.5% Methanol in water	6
5% Ethidium Bromide	6
36% Formaldehyde	3
50% Gluteraldehyde	6
0.1% Phenol	6
10-13% Sodium Hypochlorite	6

Permeation	Measured	
Performance Level	Breakthrough time (minutes)	
1	>10	
2	>30	
3	>60	
4	>120	
5	>240	
6	>480	

EN388:2003 0-0-0-0







200 ambidextrous pieces by weight













Manufactured in accordance with Quality Management System ISO 13485 & 9001
Designed and Produced to meet the applicable requirements of EN455 (Parts 1, 2, 3 and 4) / EN374 (Parts 2 & 3)



MEDIUM (7-8)

NITREX

Extra Sensitive



200 MEDIUM (7-8) MEDICAL GLOVES

GN01M

200 MEDIUM (7-8) MEDICAL GLOVES

