

GN01

NITREX Extra Sensitive



Latex-Free

Ultimate Tactile Sensitivity

Strong Chemical Resistance

Passes Viral Penetration Tests (F1671)

Tested for use with Chemotherapy Drugs

NITREX Extra Sensitive examination gloves are manufactured from a high quality nitrile formulation resulting in a strong performing glove with a high level of tactile sensitivity.

The advanced thin film technology utilised in the manufacture process results in a glove, which provides outstanding barrier protection against viruses, bacteria, chemicals and has been tested for permeation performance against 13 chemotherapy drugs.

## Specification

| Glove details | Specification              |
|---------------|----------------------------|
| Material      | Nitrile                    |
| Length        | Min. 240mm                 |
| Protein Level | Contains no latex proteins |
| Surface       | Micro-textured fingers     |
| Shape         | Ambidextrous               |
| Colour        | Blue                       |
| Sterilisation | Non sterile                |
| Shelf Life    | 3 Years                    |

## Physical Properties

| Glove details            | Properties |
|--------------------------|------------|
| Force at Break (Newtons) | >6         |
| Freedom from Holes       | AQL 1.5    |
| Length (mm)              | Min. 240   |
| Cuff Thickness (mm)      | 0.04       |
| Palm Thickness (mm)      | 0.05       |
| Finger Thickness (mm)    | 0.06       |

AQL (Acceptable Quality Level) refers to the maximum number of defective products that could be considered acceptable during the random sampling of an inspection, in this case freedom from holes in gloves. The lower the number, the fewer the holes and the higher the quality of gloves.

Sizing

| Size        | Palm width |
|-------------|------------|
| Extra Small | 75 ± 5mm   |
| Small       | 85 ± 5mm   |
| Medium      | 95 ± 5mm   |
| Large       | 106 ± 5mm  |
| Extra Large | 116 ± 5mm  |

Quantities

| Code  | Box Qty    | Case Qty |
|-------|------------|----------|
| GN01E | 200 Gloves | 10 Boxes |
| GN01S | 200 Gloves | 10 Boxes |
| GN01M | 200 Gloves | 10 Boxes |
| GN01L | 200 Gloves | 10 Boxes |
| GN01X | 200 Gloves | 10 Boxes |

Re-order Codes

| Size        | MPC   | NHS Code |
|-------------|-------|----------|
| Extra Small | GN01E | FTE1744  |
| Small       | GN01S | FTE1745  |
| Medium      | GN01M | FTE1746  |
| Large       | GN01L | FTE1747  |
| Extra Large | GN01X | FTE1748  |

**Quality standards:** Manufactured in a facility holding ISO 9001 and ISO 13485.

**Product standards:** The product is tested and manufactured to the following standards: EN455-1, EN455-2, EN455-3, EN455-4.

**Storage:** Store in a cool, dry place away from sources of heat or direct sunlight and shielded from ozone and UV light.

**Disposal:** Dispose of gloves as clinical waste. Collection case and transit case can be recycled as paper or disposed of as clinical waste. Please follow local regulation and guidance.

**Shelf life:** Three (3) years from date of manufacture.

**Country of origin:** Malaysia

**Notice:** This product is free from natural rubber latex and is suitable for use in latex-free environments. Chemical accelerators are used in the manufacture process, which can cause skin reactions in some users. Stop use of the product if skin reactions occur.



ASTM  
F1671

Find out more at [www.medicareproducts.com](http://www.medicareproducts.com)