


LATEX

Sterile
Nitrile Coating



A glove that works as hard as you do

- **Versatile glove** to be used in a wide variety of surgical environments
- **Brown-tint** blends with the wearer's skin for protection that is unobtrusive and less noticeable⁴
- Interlocking, beaded cuff design helps to **prevent roll-down**
- Proprietary hand mold with an independent thumb design allows for an **anatomical fit and natural movement**



Backed by the expertise of Cardinal Health, Protexis® Latex Surgical Gloves meet all relevant FDA and ASTM standards, including those for physical dimensions⁵, physical properties⁵ and freedom from holes⁶. Documentation is available upon request.

Protexis® Latex Surgical Gloves provide the comfort of latex with a visually appealing brown tint.



PRODUCT INFORMATION | LATEX

Properties (before aging)⁷

Tensile strength	Min 24 MPa ⁵
Stress at 500% elongation (modulus)	Max 5.5 MPa ⁵
Ultimate elongation (elasticity)	Min 750% ⁵
Puncture resistance (cuff)	Min 5N ⁸
Freedom from holes ⁶	0.65 AQL ⁵



Catalog no.	Size	Length	Thickness ⁵			Material	Color	Cuff type	Protein content	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72NS55X	5.5	11.1 in./ 282 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer		Beaded/ Rolled	Less than 50mg/dm ²	50	200
2D72NS60X	6										
2D72NS65X	6.5										
2D72NS70X	7	11.6 in./ 295 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer		Beaded/ Rolled	Less than 50mg/dm ²	50	200
2D72NS75X	7.5										
2D72NS80X	8										
2D72NS85X	8.5										
2D72NS90X	9										

1 Tested in accordance with ASTM D 5712
 2 Tested in accordance with ASTM D 6124
 3 Caution: This product contains natural rubber latex which may cause allergic reactions.
 4 As compared to other colored gloves in the Protexis® Latex Surgical Gloves Portfolio
 5 In accordance with ASTM D 3577
 6 Tested in accordance with ASTM D 5151
 7 Caution: Safe use of these gloves by or on latex-sensitized individuals has not been established.
 8 Tested in accordance with AS/NZS 4179, average test result = 8.4N
 9 93% reduction of materials used as compared to previous Cardinal Health packaging design. Cardinal Health Technical Council Conference 2011, presentation on file.

THE CARDINAL HEALTH™ PROTEXIS® SURGICAL GLOVES PROMISE:

We protect so you can perform.

As a leader in the industry with more than 50 years of surgical gloves experience, Cardinal Health is dedicated to providing protection, performance and expertise so wearers can perform confidently and focus on their patients.



Help maximize storage space: Half-fold packaging design reduces packaging material⁹



Storage recommendations: It is recommended that gloves are properly stored away from light and extreme temperatures. Gloves should also be protected from direct exposure to ozone-generated devices such as fluorescent lights, electrical motors and x-ray devices.



Expiration: 35 months from date of manufacture. Expiration date is printed on packaging.



cardinalhealth.com/surgicalgloves
 GMB-CAH-ProductSupport@cardinalhealth.com
 Customer service: 800.964.5227